





**Rocket® IPC** is indicated for intermittent, long term drainage of symptomatic, recurrent, pleural effusion, unresponsive malignant pleural effusions and abdominal ascites.






**Rocket® IPC** offers an alternative long term method of managing pleural and peritoneal ascites allowing patients to be managed in the community which offers significant advantages:

- Early discharge
- Minimised repeat admission.
- Improved quality of life
- Improved control over discomfort and breathlessness
- Reduced cost

#### Unique Flow Control:

- Gives patients or carer total control over fluid flow
- Reduces risk of accidental vacuum release

Codes	Product	Dressing Pack Contents		
<b>R54400</b> <b>FET1544</b>	<b>Rocket® IPC™ Dressing Pack with 500ml bottle</b> Supplied sterile, in cartons 5 packs	500ml pre-evacuated bottle	<ul style="list-style-type: none"> <li>• 1pr gloves (non-latex)</li> <li>• 1 x field drape</li> <li>• 4 x alcohol swabs</li> <li>• 1 x blue catheter clip - for emergency use</li> <li>• 1 x valve cap</li> <li>• 1 x foam catheter pad</li> <li>• 1 x clear-film dressing</li> <li>• 8 x gauze pads</li> <li>• 1 x Tiger Stripped waste bag</li> </ul>	
<b>R54411</b> <b>FSW706</b>	<b>Rocket® IPC™ Dressing Pack with 1000ml bottle</b> Supplied sterile, in cartons 5 packs	1000ml pre-evacuated bottle		

Code	Product Description	Code	Product Description
<b>R54410-00-CP</b> <b>FET1550</b>	<b>Rocket® IPC™ Valve Cap</b> Sterile, for single use in cartons 10 caps 	<b>R54410-00-DL</b> <b>FET1547</b>	<b>Rocket® IPC™ Drainage Line</b> Sterile for single use, in cartons 10 lines 
<b>R54411-00-00</b>	<b>Rocket® IPC™ 1000ml Drainage Bottle</b> Sterile, for single use in cartons 5 bottles 	<b>R54400-00-00</b>	<b>Rocket® IPC™ Dressing Pack</b> Sterile for single use, in cartons 5 packs 
<b>R54410-00-VV</b> <b>FSW569</b>	<b>Rocket® IPC™ Replacement Valve Kit</b> Supplied individually, sterile for single use.	 denotes 'Available on Prescription'	